



Peter C. Iwen, Ph.D., D (ABMM),  
Laboratory Director  
985900 Nebraska Medical Center  
Omaha, NE 68198  
Main: 402-559-9444  
24/7 Pager: 402-888-5588

---

# NPHL ALERT

---

Date: 12/13/2022

See below for two manufacturer's notices regarding false positive and false negative influenza results

The **Abbott Alinity** test platform has had reports of **false positive influenza B**. Abbott has resolved the issue and has contacted labs using their platform. If your lab is using the Abbott Alinity test platform for influenza testing and have not been contacted by Abbott, please contact your Abbott representative.

Roche Diagnostics has recognized **false negative influenza A** results on the **COBAS® 5800/6800/8800 Systems**. The false negative result is due to a genetic mutation in the influenza A(H1N1)pdm09 virus. Roche Diagnostics issued a global notification to affiliate organizations in November. Please contact your Roche representative with any further questions.